



Application Guide

Abri-San (Shaped) – Lying position



1. Fold the pad lengthways with the backsheet turned outwards. Do not touch the inside of the pad.

2. Pass the folded pad from front to back, creating a channel shape between the legs. If necessary, place a towel between the knees to ease the process.



3. Pull out the back of the pad horizontally, and position the upper rear at the end of the coccyx.

4. Pull up the fitting pants.



5. With the wearer lying on their back fan out the pad on the front. The pad should create a channel shape in the groin area and fit snugly to the body.



6. Adjust the fitting pants at the front. If possible, the pant legs should be placed on the backsheet edges of the pad.